RIDASCREEN[®] FAST Mandel / Almond

Art. Nr. R6901

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RIDASCREEN[®]FAST Mandel / Almond

Brief information

RIDASCREEN[®]FAST Mandel / Almond (Art. No. R6901) is a sandwich enzyme immunoassay for the quantitative analysis of almond or parts of almond in breakfast cereals, cookies, ice cream and chocolate.

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit.

The test kit is sufficient for 48 determinations (including standards).

A microtiter plate spectrophotometer is required for quantification.

Sample preparation:	homogenization, extraction and centrifugation
Time requirement:	sample preparation (for 10 samples) approx. 20 min test implementation (incubation time) 30 min
Limit of detection:	1.2 mg/kg (ppm) almond corresponding to 0.00012 % almond
Limit of quantification:	2.5 mg/kg (ppm) almond corresponding to 0.00025 % almond
Specificity:	The antibodies specifically detect proteins from almonds. Cross reaction to apricot stone > 100 % Cross reaction to mahaleb cherry > 92 % Cross reaction with other plants from genus Prunus possible. No cross-reaction with cashew nut, brazil nut, pecan, hazelnut, coconut, macadamia nut, walnut, chest nut, sunflower seeds, sesame and Lima beans.

Cross reactivities of the used antibodies have been determined for the pure food (e.g. corn flour). In a composed / processed food (e.g. corn bread) cross

reactivities might be different. Interfering substances (e.g. polyphenols) can be detected by spike experiments.

In order to increase the quality of assessment when performing ELISA-procedures, we refer additionally to our Good ELISA Practice (GEP) – Manual in the respective version. Such lists minimum standards concerning the framework conditions when using test-kits of R-Biopharm AG and performing ELISA-analysis. The manual can be retrieved, printed and downloaded under the website www.r-biopharm.com/products/food-feed-analysis.

Related products

Bioavid Lateral Flow Mandel/Almond (BL601-10 und -25)

SureFood® ALLERGEN ID Almond (S3104)

1. Intended use

RIDASCREEN®FAST Mandel / Almond is a sandwich enzyme immunoassay for the quantitative analysis of almond or parts of almond in breakfast cereals, cookies, ice cream and chocolate.

2. General

The allergen can be present as an ingredient or as a contamination in raw and cooked products. According to the **regulation (EU) No. 1169/2011**, almond must be declared on food labels as it can induce allergic reactions. Similar regulations exist e.g. in the USA, Canada, Australia and New Zealand.

Almond (*Amygdalys communis*) is one of the most popular tree nuts worldwide. Almonds are components of a lot of snacks and breakfast cereals. As ingredients almonds occur in a lot of foods especially in products of the industry for cookies and chocolate.

Almonds belong to the sub family of *Prunoideae* and are closely related to peach and apricots. The high protein content of approx. 25 % is not only nutritious but also an allergic potential. 95 % of the proteins are water soluble.

Even consumption of few milligrams of almond can induce allergic reactions in highly allergic persons. Symptoms range from allergic skin reactions up to anaphylactic shock. Since almond allergy is persistent during life and treatment of this allergy is not possible, avoiding almonds is extremely important for these patients.

3. Test principle

The basis of the test is the antigen-antibody reaction. The wells of the microtiter strips are coated with specific antibodies against almond proteins. By adding the standard or sample solution to the wells, present almond protein will bind to the specific capture antibodies and an antibody-antigen-complex is formed.

In a washing step components not bound are removed. Then antibody conjugated to peroxidase (enzyme conjugate) is added. This antibody conjugate is bound to the antibody-antigen-complex. An antibody-antigen-antibody-complex (sandwich) is formed. Any unbound conjugate is then removed in a washing step. The detection takes place by adding substrate/chromogen. The enzyme conjugate converts the chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorbance is proportional to the almond concentration in the sample.

4. Reagents provided

Each kit contains sufficient materials for 48 measurements (including standard analyses).

Component	Cap color	Format		Volume
Microtiter plate	-	Ready to use		48 wells
Allergen extraction buffer	Green	Concentrate	10x	100 ml
Standard 1*	Transparent	Ready to use	0.0 ng/ ml	1.3 ml
Standard 2*	Transparent	Ready to use	2.5 ng/ ml	1.3 ml
Standard 3*	Transparent	Ready to use	5.0 ng/ ml	1.3 ml
Standard 4*	Transparent	Ready to use	10.0 ng/ ml	1.3 ml
Standard 5*	Transparent	Ready to use	20.0 ng/ ml	1.3 ml
Wash buffer	Brown	Concentrate	10x	100 ml
Conjugate	Red	Concentrate	11x	0.7 ml
Substrate/Chromogen Red Chromogen Pro	Brown	Ready to use		10 ml
Stop solution	Yellow	Ready to use		14 ml

*) The dilution factor 20 for the sample preparation has already been considered when labeling. Therefore, the almond concentrations of samples can directly be read from the standard curve.

5. Reagents required but not provided

5.1. Equipment:

- microtiter plate spectrophotometer (450 nm)
- centrifuge + centrifugal vials
- shaker
- water bath
- laboratory mincer / grinder, pestle and mortar, Ultra-Turrax or homogenisator
- graduated pipettes
- variable 20 µl - 200 µl and 200 - 1000 µl micropipettes

5.2. Reagents:

- distilled or deionized water

6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at www.r-biopharm.com.

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The reddish substrate/chromogen is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after the expiration date on the kit label.

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- any bluish coloration of the red stained substrate/chromogen prior to test implementation

–a value of less than 0.8 absorbance units ($A_{450\text{ nm}} < 0.8$) for standard 5

9. Preparation of Samples

Almond residues of previous analyses have to be removed completely! Therefore, samples should be prepared in new or very well cleaned glass vials.

Tools such as mincers have to be cleaned thoroughly after use, in order to avoid spreading traces of almond.

The **Allergen Extraction Buffer** is provided as a 10fold concentrate. Before dilution of the buffer concentrate dissolve any crystals in a water bath at 37 °C (99 °F) completely and mix well. After that dilute the heated buffer concentrate 1:10 (1+9) with distilled water before use (i.e. 100 ml buffer concentrate + 900 ml distilled water). The diluted buffer is stable at 2 - 8 °C (35 - 46 °F) for approx. twelve weeks.

9.1. Sample preparation (**without** skim milk powder)

- 5 g of the sample should be ground well and thoroughly mixed
- weigh 1 g of sample and add 20 ml diluted Allergen Extraction Buffer (the Allergen Extraction Buffer should have already been heated to approx. 60 °C (140 °F))
- mix intensively and extract for 10 min at 60 °C (140 °F) by shaking from time to time
- centrifuge: 10 min / 2500 g / if possible at 4 °C (39 °F) and/or filter the extract (alternatively 2 ml of the extract can be centrifuged with high speed for 10 min in reaction caps by using a microcentrifuge)
- use 100 µl of the supernatant or filtrate per well in the assay

9.2. Sample preparation (**with** skim milk powder) for spices and cacao

- 5 g of the sample should be ground well and thoroughly mixed
- weigh 1 g of sample and add 1 g skim milk powder
- add 20 ml diluted Allergen Extraction buffer (the Allergen Extraction buffer should have already been heated to approx. 60 °C (140 °F))
- mix intensively and extract for 10 min at 60 °C (140 °F) by shaking casually
- centrifuge: 10 min / 2500 g / if possible at 4 °C (39 °F) and/or filter the extract (alternatively 2 ml of the extract can be centrifuged with high speed for 10 min in reaction caps by using a microcentrifuge)
- use 100 µl of the supernatant or filtrate per well in the assay

Remark:

The sample extracts can be stored at 2 - 8 °C (35 - 46 °F) for 3 days.

10. Test implementation

10.1. Preliminary comments

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **Conjugate** (bottle with red cap) is provided as a concentrate. Since the diluted Conjugate has a limited stability, only the amount which actually is needed should be reconstituted. Before pipetting, the Conjugate should be shaken carefully. For reconstitution, the Conjugate is diluted 1:11 (1+10) in distilled water (e.g. 200 µl Conjugate + 2 ml distilled water, sufficient for 2 microtiter strips).

The **Washing Buffer** is provided as a 10fold concentrate. Prior to dilution, dissolve eventually formed crystals by incubating the buffer in a water bath at 37 °C (99 °F). Before use, the buffer has to be diluted 1:10 (1+9) with distilled water (i.e. 100 ml buffer concentrate + 900 ml dist. water). The diluted buffer is stable at 2 - 8 °C (35 - 46 °F) for approx. four weeks.

10.2. Test procedure

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

Do not use more than 3 strips (24 wells) at a time. In the case of more than three strips, a second uncoated plate (e.g. low binding from Greiner bio-one Cat.-No. 655101) should be used as a pre-plate to avoid a time shift over the microtiter plate. All standards and samples are pipetted into the uncoated plate (at least 150 µl per well) and then quickly transferred to the coated microtiter plate with an 8-channel pipette.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run. Record standard and sample positions.
2. Add 100 µl of each standard or prepared sample to separate wells and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).
3. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl washing buffer (see 10.1.) and pour out the liquid again. Repeat two more times.

4. Add 100 µl of the diluted enzyme conjugate (see 10.1.) to each well. Mix gently by shaking the plate manually and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).
5. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl washing buffer (see 10.1.) and pour out the liquid again. Repeat two more times.
6. Add 100 µl of the reddish substrate/chromogen to each well. Mix gently by shaking the plate manually and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
7. Add 100 µl of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 10 minutes after addition of stop solution.

11. Results

A special software, the RIDA[®]SOFT Win / RIDA[®]SOFT Win.net (Art. Nr. Z9996), is available for evaluation of the RIDASCREEN[®] enzyme immunoassays. The calculation should be done by use of a cubic spline function. The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit.

In comparison with the certificate, higher values of the absorbance (A450 nm) for the calibration curve, especially for the zero standard, may be a result of insufficient washing or allergen contamination.

A further dilution and new detection of the samples is recommended for absorbance values (A450 nm) > standard 5.

Please note:

When working in accordance with the regulation stated, the dilution factor is 20. The allergen concentration can be read directly from the standard curve (see 4. *) - the sample dilution factor of 20 is already taken into account.

For sample dilutions of more than 1:20 the further dilution factor must be considered for the calculation of the almond concentration.

In general:

Samples tested negative still could contain an allergen contamination below the limit of detection of the assay, or they might contain other allergen components like lipids for example.

Due to the multitude of food types, matrix effects cannot be excluded. In processed food, proteins may be altered or fragmented, this may have an impact on the recovery.

Recommendation:

In order to ensure a high analytical performance:

- Each sample material should be analyzed in duplicates
- Use also allergen-free and allergen- containing (spiked) samples as test controls
- Carry out spiking experiments because due to the multitude of food types, matrix effects cannot be excluded.
- For details using the ChemWell® or GEMINI automation please contact sales@r-biopharm.de.

The product information folder with further information is available on request from your local distributor or R-Biopharm AG.

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG